

Elite™ Hydrogen Peroxide Detection Kit (Red Fluorescence)

CATALOG NUMBER: CA-H501, 500 assays

Description

Hydrogen peroxide (H_2O_2) is a reactive oxygen metabolic by-product that serves as a key regulator for a number of oxidative stress-related states. It is involved in several biological events that have been linked to asthma, atherosclerosis, diabetic vasculopathy, osteoporosis, a number of neurodegenerative diseases and Down's syndrome. Perhaps the most intriguing aspect of H_2O_2 biology is the recent report that antibodies have the capacity to convert molecular oxygen into hydrogen peroxide to contribute to the normal recognition and destruction processes of the immune system. Measurement of this reactive species will help to determine how oxidative stress modulates a variety of intracellular pathways.

The Elite™ Fluorimetric Hydrogen Peroxide Detection Kit uses our non-fluorescent Elite™ Red peroxidase substrate to quantify hydrogen peroxide in solutions and cell extracts. It can also be used to detect a variety of oxidase activities through enzyme-coupled reactions. The kit is an optimized “mix and read” assay that is compatible with HTS liquid handling instruments. It provides a sensitive, one-step fluorometric assay to detect as little as 3 picomoles of H_2O_2 in a 100 μ L assay volume (30 nM, Figure 1). The assay can be performed in a convenient 96-well or 384-well microtiter-plate format and readily adapted to automation. Its signal can be easily read by either a fluorescence microplate reader at Ex/Em = ~540/590 nm or an absorbance microplate reader at ~570 nm.

Kit Components

• Component A: Elite™ Red Peroxidase Substrate	1 vial
• Component B: H_2O_2	1 vial (3% stabilized solution, 200 μ l)
• Component C: Assay Buffer	1 bottle (100 ml)
• Component D: Horseradish Peroxidase	1 vial (20 units)
• Component E: DMSO	1 vial (1 ml)

Storage

Keep in freezer (-20 °C) and avoid exposure to light. All Components are stable for at least 6 months after receipt if stored properly.

Materials Required (but not supplied)

- 96 microplates: Tissue culture microplate with black wall and clear bottom is recommended.
- Fluorescence microplate reader.

Assay Protocol (for 96-well plate)

1. Preparation of stock solutions:

- 1.1. 100X Elite™ Red peroxidase substrate stock solution: Add 250 μ L of DMSO (**Component E**) into the vial of Elite™ Red Substrate (**Component A**). The stock solution should be used promptly; any remaining solution should be aliquoted and refrozen at -20 °C.

Note: Avoid repeated freeze-thaw cycles and protect from light.

- 1.2. 20 U/mL Peroxidase stock solution: Add 1 mL of Assay Buffer (**Component C**) into the vial of Horseradish Peroxidase (**Component D**).

Note: The unused HRP solution should be divided into single use aliquots and stored at -20 °C.

- 1.3. 20 mM H_2O_2 stock solution: Add 22.7 μ L of 3% H_2O_2 (0.88 M, Component B) into 977 μ L of Assay Buffer (**Component C**).

Note: The diluted H_2O_2 solution is not stable. The unused portion should be discarded.



2. Preparation of H₂O₂ reaction mixture:

Prepare the H₂O₂ reaction mixture per table 1 and keep from light:

Table 1. H₂O₂ Reaction mixture for one 96-well plate (2X)

Components	Volume
Elite™ Red Peroxidase Substrate Stock Solution (100x, from Step 1.1)	50 µl
20 U/ml Peroxidase Stock Solution (from Step 1.2)	200 µl
Assay Buffer (Component C)	4.75 ml
Total Volume	5 ml

3. Preparation of serial dilutions of H₂O₂ standard (0 to 10 µM):

Warning:

- **The component A** is unstable in the presence of thiols such as DTT and β-ercaptoethanol. Thiols higher than 10 µM (final concentration) would significantly decrease the assay dynamic range.
- NADH and glutathione (reduced form: GSH) may interfere with the assay.

3.1. Add 1 µL of 20 mM H₂O₂ solution (from Step 1.3) into 1999 µL of Assay Buffer (**Component C**) to get a 10 µM H₂O₂ standard.

3.2. Take 200 µL of 10 µM H₂O₂ standard to perform 1:3 serial dilutions to get 3, 1, 0.3, 0.1, 0.03, 0.01 and 0 µM serial dilutions of H₂O₂ standard.

3.3. Add serial dilutions of H₂O₂ standard and H₂O₂-containing test samples into a solid black 96-well microplate as described in Tables 2 and 3.

Table 1. Layout of H₂O₂ standards and test samples in a solid black 96-well microplate

BL	BL	TS	TS						
HS1	HS1								
HS2	HS2										
HS3	HS3										
HS4	HS4										
HS5	HS5										
HS6	HS6										
HS7	HS7										

Note: HS= H₂O₂ Standards; BL=Blank Control; TS=Test Samples.

Table 2. Reagent composition for each well

H ₂ O ₂ Standards	Blank Control	Test Sample
Serial Dilutions: 50 µl	Assay Buffer (Component C): 50 µl	50 µl

Note: Add the serially diluted H₂O₂ standards from 0.01 µM to 10 µM into wells from HS1 to HS7 in duplicate. High concentration of H₂O₂ (e.g., > 100 µM, final concentration) may cause reduced fluorescence signal due to the overoxidation of Elite™ Red (to a non-fluorescent product).



4. Run H₂O₂ assay in supernatants reaction:

- 4.1. Add 50 µL of H₂O₂ reaction mixture (from Step 2) into each well of H₂O₂ standard, blank control, and test samples (see Step 3.3) to make the total H₂O₂ assay volume of 100 µL/well.

Note: For a 384-well plate, add 25 µL of sample and 25 µL of H₂O₂ reaction mixture in each well. Incubate the reaction mixture at room temperature for 15 to 30 minutes (protected from light).

- 4.2. Incubate the reaction at room temperature for 15 to 30 minutes, protected from light.

- 4.3. Monitor the fluorescence increase with a fluorescence plate reader at Ex/Em = 540 ± 10 / 590 ± 10 nm (optimal Ex/Em = 540/590 nm).

Note: The contents of the plate can also be transferred to a white clear bottom plate and read by an absorbance microplate reader at the wavelength of 576 ± 5 nm. The absorption detection has lower sensitivity compared to fluorescence reading.

5. Run H₂O₂ assay for cells:

The Elite™ Hydrogen Peroxide Assay Kit can be used to measure the release of H₂O₂ from cells. The following is a suggested protocol that can be modified to meet the specific research needs.

- 5.1. The H₂O₂ reaction mixture should be prepared as Step 2 except that the Assay Buffer (**Component C**) should be replaced with the media that is used in your cell culture system. Suggested media including (a) Krebs Ringers Phosphate Buffer (KRPB); (b). Hanks Balanced Salt Solution (HBSS); or (c) Serum-free media.

- 5.2. Prepare cells in a 96-well plate (50 - 100 µL/well), and activate the cells as desired.

Note: The negative controls (media alone and non-activated cells) are included for measuring background fluorescence.

- 5.3. Add 50 µL of H₂O₂ reaction mixture (from Step 5.1) into each well of cells and H₂O₂ standards (from Step 3.3).

Note: For a 384-well plate, add 25 µL of cells and 25 µL of H₂O₂ reaction mixture into each well.

- 5.4. Incubate the reaction at room temperature for 15 to 30 minutes, protected from light.

- 5.5. Monitor the fluorescence increase with a fluorescence plate reader at Ex/Em = 540 ± 10 / 590 ± 10 nm (optimal Ex/Em = 540/590 nm).



Data Analysis

The fluorescence in blank wells (with the assay buffer only) is used as a control, and is subtracted from the values for those wells with the H₂O₂ reactions. A H₂O₂ standard curve is shown in Figure 1.

Note: The fluorescence background increases with time, thus it is important to subtract the fluorescence intensity value of the blank wells for each data point.




Figure 1. H₂O₂ dose response was measured in a 384-well black plate with the Elite™ Hydrogen Peroxide Detection Kit using a Gemini fluorescence microplate reader (Molecular Devices). As low as 0.03 μM H₂O₂ can be detected with 30 minutes incubation (n=3). The insert shows the low levels of H₂O₂ detection.

References

- Hoffmann, O. M., Becker, D., and Weber, J. R. (2007) J Cereb Blood Flow Metab.
- Funk, R. S., and Krise, J. P. (2007) Mol Pharm 4, 154-9.
- Krebs, B., Wiebelitz, A., Balitzki-Korte, B., Vassallo, N., Paluch, S., Mitteregger, G., Onodera, T., Kretzschmar, H. A., and Herms, J. (2007) J Neurochem 100, 358-67.
- Yang, Y., Xu, S., An, L., and Chen, N. (2007) J Plant Physiol.
- Lee, J. E., Kim, H., Jang, H., Cho, E. J., and Youn, H. D. (2007) J Neurochem.
- Goth, L. (2006) Redox Rep 11, 281-2.
- Laloi, C., Stachowiak, M., Pers-Kamczyc, E., Warzych, E., Murgia, I., and Apel, K. (2007) Proc Natl Acad Sci U S A 104, 672-7.