

Elite™ High Fidelity DNA Polymerase

Catalog No. DP-003-0250, 250 Units

DP-003-0500, 500 Units

APPLICATION	Due to its 3'-5' proofreading exonuclease activity, Elite™ high fidelity DNA polymerase is strongly suggested for extremely high fidelity PCR applications, such as subcloning of PCR amplified cDNA, cDNA library construction, etc.
DESCRIPTION	Elite™ High Fidelity DNA Polymerase is a thermostable enzyme possessing 5'-3' DNA polymerase and 3'-5' proofreading exonuclease activities. It is isolated from the hyperthermophilic marine archae <i>Pyrococcus furiosus</i> (Pfu). Base misinsertions that may occur during polymerization are rapidly excised by the proofreading activity of the polymerase. Elite™ DNA Polymerase generates blunt-ended PCR fragments. Elite™ DNA Polymerase can amplify fragments up to 5 kb at extremely high fidelity. A mixture of BioTherm™ DNA Polymerase with Elite™ DNA Polymerase provides more robust synthesis of longer (>5 kb) amplification products.
CONCENTRATION	5 units/μl.
UNIT DEFINITION	One unit is defined as the amount of enzyme that incorporates 10 nmoles of dNTPs into acid-insoluble form in 30 minutes at 72°C under the assay conditions (25 mM TAPS (tris-(hydroxymethyl)methyl-amino-propane-sulphonic acid, sodium salt) pH 9.3 (at 25°C), 50 mM KCl, 2 mM MgCl ₂ , 1 mM β-mercaptoethanol) and activated calf thymus DNA as substrate.
STORAGE BUFFER	10 mM K-phosphate buffer pH 7.0, 100 mM NaCl, 0.5 mM EDTA, 1 mM DTT, 0.01% Tween 20; 50% glycerol (v/v)
5X REACTION BUFFER	83 mM (NH ₄) ₂ SO ₄ , 350 mM Tris-HCl pH 8.8 (at 25°C), 22 mM MgCl ₂ , 0,75% Triton X100, 100 mM KCl.
STORAGE TEMPERATURE	Store at -20°C in a constant temperature freezer.

REFERENCES

1. Fiedler, et al. (2008). *J Bacteriol* 190: 6559—6567.
2. Jahromi, et al. (2009). *Hybridoma (Larchmt)*. 28: 305--313.
3. Moreno, et al. (2008). *J Biol Chem*. 283: 9633--9641.



GUIDELINE OF PCR AMPLIFICATION:

Set-up of a Standard PCR Amplification	
Component	Amount
ddH ₂ O	- μl
10x buffer	10 μl
dNTP (25 mM)	0.8 μl
DNA template	1 ng-1 μg
Primer 1 (20 μM)	5 μl
Primer 2 (20 μM)	5 μl
Taq DNA polymerase (5 units/μl)	0.4-2 μl
Total volume	100 μl

Note:

1. Mg²⁺ concentration varies in different PCR reaction buffer system. For amplification of genomic DNA sequence, the optimal Mg²⁺ should be determined by adding Mg²⁺ to a final concentration of 1.5 mM to 4.5 mM;
2. The amount of DNA template required varies depending on the type of DNA being amplified. Generally, 50 ng to 1000 ng of genomic DNA is recommended; Less DNA (1 ng-100 ng) can be used for amplification of plasmid DNA, purified DNA, and virus DNA;
3. Primer concentration between 0.2 μM and 1 μM are recommended (approximately 100 ng to 250 ng for typical 18- to 25-mer oligonucleotide primers in 100 μl reaction volume).
4. The amount of Taq DNA polymerase varies depending on the length of template to be amplified. Successful amplification can usually be achieved using 2-5 units of enzyme/100 μl reaction volume (higher for long template amplification).

Standard PCR Temperature cycling Program		
1. 94°C	4 min	
2. 94°C	30-60 sec	25-35 Cycles
55-65°C	30-60 sec	
72°C	1 min (add 1 min per kb target sequence length)	
3. 72°C	5 min	
4. 4°C		

Note:

1. High quality thin-wall PCR tubes, strips, and plates are highly recommended for getting a consistent and duplicable PCR results. This is particularly important when amplification volume is less than 20 μl. EU thin-wall tubes can produce the most reliable results and the amplification volume can be as low as 5 μl.
2. The typical annealing temperature is between 55 °C and 72 °C. The optimal annealing temperature could be calculated by "T_m - 5 °C".