

EU Plastics PCR Compatibility Chart



For Real-Time Cyclers see page 17			PCR Tubes						PCR Strips						PCR Plates																																			
Manufacturer	Model	Wells	C(F)78401	C(F)79401	B79401	B79001	FA78801	B/C79801	C(F(A))79601	C(F)78601	B79901	B79601	B72711	C(F)78201	B79201	B76601	B70501	AB70651	B60101	B70651	B70671	AB17500	AB19800	B71501																										
ABI	GeneAmp [®] PCR System 2400	24	☑	☑	☑	☑	☑	☑	☑	☑	☑	☐	☐	☑	☑	☐	✗	☑	☑	☑	☑	☑	☑	☑																										
ABI	GeneAmp [®] PCR System 9600	96	☑	☑	☑	☑	☑	☑	☑	☑	☑	☐	☐	☑	☑	☐	☑	☐	☐	☑	☑	☑	☑	☑																										
ABI	2700 Thermal Cycler	96	☑	☑	☑	■	☑	☑	☑	☑	☑	☐	☐	☑	☑	☐	☑	☐	☐	☑	☑	☑	☑	☑																										
ABI	2720 Thermal Cycler	96	☑	☑	☑	■	☑	☑	☑	☑	☑	☐	☐	☑	☑	☐	☑	☐	☐	☑	☑	☑	☑	☑																										
ABI	GeneAmp [®] PCR System 9700	60	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑																										
ABI	GeneAmp [®] PCR System 9700	96	☑	☑	☑	■	☑	☑	☑	☑	☑	☐	☐	☑	☑	☐	☑	☐	☐	☑	☑	☑	☑	☑																										
ABI	GeneAmp [®] PCR System 9700	384	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑																										
ABI	9800 Fast Thermal Cycler	96 fast	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑																										
Manufacturer			Model																								Wells																							
Bioer	XP Cycler	60	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑																										
Bioer	XP Cycler	30/48/96/384	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑																										
Bioer	Little Genius	16	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑																										
Bioer	Little Genius	25	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑																										
Bioer	Life Express	48 (0.5ml)	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑																										
Bioer	Life Express	48 (0.2ml)	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑																										
Bioer	Life Express	96	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑																										
Manufacturer			Model																								Wells																							
Biometra	Uno	24/48/96	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑																										
Biometra	Uno II	24/48/96/384	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑																										
Biometra	Tpersonal	20/48	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑																										
Biometra	T1 (plus)	48/96/384	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑																										
Biometra	TGradient	48/96/384	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑																										
Biometra	T3/T3000	20/48	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑																										
Biometra	TProfessional	96/384	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑																										
Biometra	TRobot	96/384	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑																										
Manufacturer			Model																								Wells																							
Bio-Rad	Gene Cycler	24	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑																										
Bio-Rad	iCycler [®]	48/96/384	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑																										
Bio-Rad	MyCycler	96	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑																										
For MJ-Research Systems see Page 15																																																		
Manufacturer			Model																								Wells																							
CLP	ATC 201	48	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑																										
CLP	ATC 201	96	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑																										
CLP	ATC 401	48	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑																										
CLP	ATC 401	96	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑																										
Manufacturer			Model																								Wells																							
Corbett Research	Palm-Cycler [®]	60/96/384	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑																										

- ☑ Fits this cycler
- ☐ Does not fit this cycler
- ☐ Requires EU adaptor to fit cycler see page 34
- Minimum of 20 tubes per block
- ✓ Most optimal plastic for this thermocycler
- ✗ Can be cut to fit this cycler

EU Plastics PCR Compatibility Chart



For Real-Time Cyclers see page 17			PCR Tubes						PCR Strips						PCR Plates																																	
Manufacturer	Model	Wells	C(F)78401	C(F)79401	B79401	B79001	FA78801	B/C79801	C(F(A))79601	C(F)78601	B79901	B79601	B72711	C(F)78201	B79201	B76601	B70501	AB70651	B60101	B70651	B70671	AB17500	AB19800	B71501																								
Eppendorf	Mastercycler	77	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒																								
Eppendorf	Mastercycler	96	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☐	☐	☒																							
Eppendorf	Mastercycler	384	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒																								
Eppendorf	Mastercycler gradient	77	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒																								
Eppendorf	Mastercycler gradient	96	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☐	☐	☒																							
Eppendorf	Mastercycler gradient	384	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒																								
Eppendorf	Mastercycler personal	16	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒																								
Eppendorf	Mastercycler personal	25	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒																								
Eppendorf	Mastercycler ep gradient	96	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☐	☐	☒																							
Eppendorf	Mastercycler ep gradient S	96	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☐	☐	☒																							
Eppendorf	Mastercycler ep 384	384	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒																								
Manufacturer			Model																							Wells																						
GRI	G-storm GSI	48 0.2ml/0.5ml	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☐	☐	☒																						
GRI	G-storm GSI	96/384	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒																						
GRI	G-storm GSX	48 0.2ml/0.5ml	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☐	☐	☒																						
GRI	G-storm GSX	96/384	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒																						
GRI	G-storm GS4	48 0.2ml/0.5ml	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☐	☐	☒																						
GRI	G-storm GS4	96/384	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒																						
Manufacturer			Model																							Wells																						
MIJ Research/Bio-Rad	PTC-150	16/25	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒																							
MIJ Research/Bio-Rad	Mini	48	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒																							
MIJ Research/Bio-Rad	PTC-100	60/96	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☐	☐	☒																						
MIJ Research/Bio-Rad	PTC-200	30/48/96/384	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☐	☐	☒																						
MIJ Research/Bio-Rad	PTC-221	30/48/96/384	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☐	☐	☒																						
MIJ Research/Bio-Rad	PTC-220	30/48/96/384	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☐	☐	☒																						
MIJ Research/Bio-Rad	PTC-225	30/48/96/384	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☐	☐	☒																						
For Bio-Rad System see Page 14																																																
Manufacturer			Model																							Wells																						
MMG	Primus 25	13/25	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒																						
MMG	Primus 96	96	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☐	☐	☒																						
MMG	Primus 384	384	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒																						
MMG	Primus HT Dualblock	96/384	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☐	☐	☒																						
MMG	Primus Multiblock System	96/384	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☐	☐	☒																						
Manufacturer			Model																							Wells																						
Stratagene	Robo-Cycler	96	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☐	☐	☒																						
Manufacturer			Model																							Wells																						
Takara	TP 100	24	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒																						
Takara	TP 240	24	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒																						
Takara	TP 600	96	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒																						
Takara	TP 650	96	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒																						
Takara	TP 3000	96	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒																						

- ☒ Fits this cycler
- ☒ Does not fit this cycler
- ☐ Requires EU adaptor to fit cycler see page 34
- Minimum of 20 tubes per block
- ✓ Most optimal plastic for this thermocycler
- ≧ Can be cut to fit this cycler

EU Plastics PCR Compatibility Chart



For Real-Time Cyclers see page 17			PCR Tubes						PCR Strips						PCR Plates										
Manufacturer	Model	Wells	C(F)78401	C(F)79401	B79401	B79001	FA78801	B/C79801	C(F(A))79601	C(F)78601	B79901	B79601	B72711	C(F)78201	B79201	B76601	B70501	AB70651	B60101	B70651	B70671	AB17500	AB19800	B71501	
Techne	Genius	40/96/384	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	□	□	✓
Techne	TouchGene	40/96	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	□	□	✓
Techne	FlexiGene/FlexiQuad	96/384	✓	✓	✓	✓	✓	✗	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	□	□	✓
Techne	Progene	20/25	✓	✓	✓	✓	✓	✓	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗
Techne	TechGene	25/96	✓	✓	✓	✓	✓	✗	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	□	□	✓
Techne	Gene-E	96	✓	✓	✓	✓	✓	✗	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	□	□	✓
Techne	TC-312	20/25	✓	✓	✓	✓	✓	✓	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗
Techne	TC-412	96	✓	✓	✓	✓	✓	✗	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	□	□	✓
Techne	TC-512	96/384	✓	✓	✓	✓	✓	✗	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	□	□	✓
Manufacturer	Model	Wells																							
ThermoHybaid	PfE	96	✓	✓	✓	✓	✓	✗	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	□	□	✓
ThermoHybaid	OmniGene	96	✓	✓	✓	✓	✓	✗	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	□	□	✓
ThermoHybaid	PCR Express	48/96/384	✓	✓	✓	✓	✓	✗	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	□	□	✓
ThermoHybaid	Omni-E	96	✓	✓	✓	✓	✓	✗	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	□	□	✓
ThermoHybaid	TouchDown	48/96	✓	✓	✓	✓	✓	✗	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	□	□	✓
ThermoHybaid	Pf2	96/384	✓	✓	✓	✓	✓	✗	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	□	□	✓
ThermoHybaid	PCR sprint	24	✓	✓	✓	✓	✓	✗	✓	✓	✓	✓	✓	✓	✓	✓	✗	✗	✗	✗	✗	✗	✗	✗	✗
ThermoHybaid	MBS	96/384	✓	✓	✓	✓	✓	✗	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	□	□	✓
ThermoHybaid	MBSR	96/384	✓	✓	✗	✗	✓	✗	✓	✓	✓	✗	✗	✓	✓	✗	✓	✓	✓	✓	✓	✓	□	□	✓

EU Plastics DNA Sequencer Compatibility Chart

Manufacturer	Model	PCR Plates											
		B70501	AB70651	B60101	B70651	B70671	AB17500	AB19800	B71501				
ABI	310	✓	✗	✗	✗	✗	✗	✗	✗				
ABI	3100	✓	✗	✗	✗	✗	✗	✗	✗				
ABI	3130	✓	✗	✗	✗	✗	✗	✗	✗				
ABI	3700	✓	✗	✗	✗	✗	✓	✗	✗				
ABI	3730	✓	✗	✗	✗	✗	✓	✗	✗				
Manufacturer	Model												
Amersham	MegaBACE™ 500	✗	✗	✗	✗	✓	✗	✗	✗				
Amersham	MegaBACE™ 1000	✗	✗	✗	✗	✓	✗	✗	✗				
Amersham	MegaBACE™ 4000	✗	✗	✗	✗	✗	✗	✗	✗				
Manufacturer	Model												
Beckman	CEQ 8000	✗	✗	✗	✗	✗	✓	✗	✗				
Beckman	CEQ 8800	✗	✗	✗	✗	✗	✓	✗	✗				
Manufacturer	Model												
Transgenomics	Wave MD System 4000	✓	✗	✗	✗	✓	✗	✗	✗				

- ✓ Fits this cycler
- ✗ Does not fit this cycler
- Requires EU adaptor to fit cycler see page 34
- Minimum of 20 tubes per block
- ✓ Most optimal plastic for this thermocycler
- ✗ Can be cut to fit this cycler

EU Plastics Real-Time PCR Compatibility Chart



For regular Cyclers see pages 14 through 16			PCR Tubes						PCR Strips						PCR Plates																																			
Manufacturer	Model	Wells	C(F)78401	C(F)79401	B79401	B79001	FA78801	B/C79801	C(F(A))79601	C(F)78601	B79901	B79601	B72711	C(F)78201	B79201	B76601	B70501	AB70651	B60101	B70651	B70671	AB17500	AB19800	B71501																										
ABI	GeneAmp [®] SDS 5700	96	✓	✓	✗	✓	✓	✗	✓	✓	✓	✗	✗	✓	✓	✗	✓	✗	✗	✓	✗	✓	✗	✓																										
ABI	Prism [®] SDS 7000	96	✓	✓	✗	✓	✓	✗	✓	✓	✓	✗	✗	✓	✓	✗	✓	✗	✗	✓	✗	✓	✗	✓																										
ABI	Prism [®] SDS 7300	96	✓	✓	✗	✓	✓	✗	✓	✓	✓	✗	✗	✓	✓	✗	✓	✗	✗	✓	✗	✓	✗	✓																										
ABI	Prism [®] SDS 7500	96	✓	✓	✗	✓	✓	✗	✓	✓	✓	✗	✗	✓	✓	✗	✓	✗	✗	✓	✗	✓	✗	✓																										
ABI	Prism [®] SDS 7500	96 fast	✗	✗	✗	✗	✗	✗	✗	✗	✗	✓	✓	✗	✗	✓	✗	✗	✗	✗	✗	✗	✓	✗																										
ABI	Prism [®] SDS 7700	96	✓	✓	✗	✓	✓	✗	✓	✓	✓	✗	✗	✓	✓	✗	✓	✗	✗	✓	✗	✓	✗	✓																										
ABI	Prism [®] SDS 7900	96	✗	✗	✗	✗	✓	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✓	✗	✗																										
ABI	Prism [®] SDS 7900	384	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗																										
ABI	Prism [®] SDS 7900	96 fast	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✓	✗																										
Manufacturer			Model																								Wells																							
Bioer	Line-Gene	33	✓	✓	✓	✓	✓	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗																									
Bioer	Line-Gene II	66	✓	✓	✓	✓	✓	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗																									
Manufacturer			Model																								Wells																							
Bio-Rad	iCycler iQ	96	✓	✓	✓	✓	✓	✗	✓	✓	✓	✗	✗	✓	✓	✗	✓	✗	✗	✓	✗	✗	✗	✓																										
Bio-Rad	iQ5	96	✓	✓	✓	✓	✓	✗	✓	✓	✓	✗	✗	✓	✓	✗	✓	✗	✗	✓	✗	✗	✗	✓																										
Bio-Rad	MyiQ	96	✓	✓	✓	✓	✓	✗	✓	✓	✓	✗	✗	✓	✓	✗	✓	✗	✗	✓	✗	✗	✗	✓																										
Manufacturer			Model																								Wells																							
Corbett Research	Rotor-Gene [®] 2000	36/72	✓	✓	✓	✓	✓	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗																										
Corbett Research	Rotor-Gene [®] 3000	36/72	✓	✓	✓	✓	✓	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗																										
Corbett Research	Rotor-Gene [®] 6000	36/72	✓	✓	✓	✓	✓	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗																										
Manufacturer			Model																								Wells																							
Eppendorf	MasterCycler [®] ep realplex	96	✓	✓	✓	✓	✓	✗	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✗	✗																										
Manufacturer			Model																								Wells																							
MJ Research/Bio-Rad	Opticon [®]	96	✗	✗	✗	✗	✗	✗	✗	✗	✗	✓	✓	✗	✗	✓	✗	✗	✓	✓	✗	✓	✗	✗																										
MJ Research/Bio-Rad	Opticon2	96	✗	✗	✗	✗	✗	✗	✗	✗	✗	✓	✓	✗	✗	✓	✗	✗	✓	✓	✗	✓	✗	✗																										
MJ Research/Bio-Rad	MiniOpticon [®]	48	✓	✓	✓	✓	✓	✗	✓	✓	✓	✓	✓	✓	✓	✓	✗	✗	✗	✗	✗	✗	✗	✓																										
MJ Research/Bio-Rad	Chromo4	96	✓	✓	✓	✓	✓	✗	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✗	✗																										
Manufacturer			Model																								Wells																							
Roche	LightCycler [®] 480	96	✓	✓	✗	✗	✓	✗	✗	✗	✗	✓	✓	✗	✗	✓	✗	✗	✓	✓	✗	✓	✗	✗																										
Roche	LightCycler [®] 480	384	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗																										
Roche	LightTyper System	96	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✓	✓	✗	✓	✗	✗																										
Roche	LightTyper System	384	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗																										
Manufacturer			Model																								Wells																							
Stratagene	Mx3000P [®]	96	✓	✓	✓	✓	✓	✗	✓	✓	✗	✗	✗	✗	✗	✗	✓	✗	✗	✗	✗	✗	✗	✓																										
Stratagene	Mx3005P [®]	96	✓	✓	✓	✓	✓	✗	✓	✓	✗	✗	✗	✗	✗	✗	✓	✗	✗	✗	✗	✗	✗	✓																										
Stratagene	Mx4000	96	✓	✓	✓	✓	✓	✗	✓	✓	✗	✗	✗	✗	✗	✗	✓	✗	✗	✗	✗	✗	✗	✓																										
Manufacturer			Model																								Wells																							
Techne	Quantica [®]	96	✗	✗	✗	✗	✗	✗	✗	✗	✗	✓	✓	✗	✗	✓	✗	✗	✓	✓	✗	✓	✗	✗																										

- ✓ Fits this cycler
- ✗ Does not fit this cycler
- ✗ Requires EU adaptor to fit cycler see page 34
- Minimum of 20 tubes per block
- ✓ Most optimal plastic for this thermocycler
- ✗ Can be cut to fit this cycler