

Phosphodiesterase 1B (PDE1B) Stable Cell Line


CATALOG NUMBER: CL-02-PDE1B

Introduction

The protein encoded by PDE1B gene belongs to the cyclic nucleotide phosphodiesterase (PDE) family, and PDE1 subfamily. Members of the PDE1 family are calmodulin-dependent PDEs that are stimulated by a calcium-calmodulin complex. This PDE has dual-specificity for the second messengers, cAMP and cGMP, with a preference for cGMP as a substrate. cAMP and cGMP function as key regulators of many important physiological processes.

Description

Human PDE1B ACTOne™ is a CHO-CNG-sGC cell line that expresses human PDE1B. CHO-CNG-sGC cells express a modified CNG (Cyclic Nucleotide Gated) channel that opens in response to elevated intracellular cAMP levels and consequently result in ion flux (often detectable by calcium-responsive dye, [Cat# CA-C155](#)) and cell membrane depolarization which can be easily measured with fluorescent Membrane Potential Dye ([Cat# CA-M165](#)). The assay allows both end-point and kinetic measurement of intracellular cAMP changes with a FLIPR, or a fluorescence microplate reader.



Parental Cells

CHO-CNG-sGC cells (originally developed by BD Biosciences) (Cat# CL-02-PC50)

Gene/Enzyme Introduced

PDE1B (Genebank No. NP_000915.1)

Applications

- cGMP dependent human PDE1B cell based assay
- cell based high-throughput screening of endogenous human PDE1B inhibitors

Functional Test

- this cell line has been tested positive for PDE1B specific response
- surviving rate: More than 2.5 million/vial on the second day after thawing
- the receptor specific activity is stable for 10 weeks continuous passage



Mycoplasma Contamination Test

This lot of cells has been tested and found to be free of mycoplasma contamination.

Content

- Stable cells (CL-02-PDE1B): 1 mL (1×10^6 cells/mL in 70% DMEM, 20% FBS, 10% DMSO)

Growth Properties

Adherent

Cell Culture Medium

- DMEM/Nutrient F-12 Ham-10% FBS supplemented with 250 μ g/ml G418, 1 μ g/ml Puromycin and 5 μ g/ml blasticidin.
- Freezing medium: 10% DMSO, 90% complete cell culture medium

Storage

Remove the frozen cells from the dry ice packaging and immediately place the cells at a temperature below -130°C , preferably in liquid nitrogen vapor, until ready for use.

Data Analysis




Figure 1. Response of ACTOne™ cGMP-PDE1B cell line & parental cell line to Vardenafil.

ACTOne™ cGMP-PDE1B cells and parental cells (Cat# CL-02-PC50) were plated overnight in 20 μ l culture medium on a 384 well microplate. The next day, cells were dye-loaded with 20 μ l/well of membrane potential dye (Cat# CA-M165). After 2 hours of incubation at room temperature, baseline was recorded using a FlexStation (Molecular Devices) (F_0). 10 μ l of PDE inhibitors at various concentrations (with 125 μM Bay41-2272) were added to the cell plate, and the data was recorded 40 minutes (F_t) after drug addition. Dose response curves were generated by Prism.

- Dose response curve of Vardenafil in ACTOne cGMP-PDE1B cell line. $\text{EC}_{50} = 252 \text{ nM}$ in the presence of 25 μM of Bay 41-2272
- Parental cells do not respond to Vardenafil in the presence of 25 μM of Bay 41-2272

Notice to Purchaser

This cell line is sold with a use license. It may not be transferred to third parties, resold, modified for resale without written approval of eEnzyme LLC. Refer to the license agreements for details on the usage restrictions.

